



Auxiliary Lens 0.5X

Item Type: Auxiliary Objective Lens

Model: Aux 0.5

Material: High-quality Aluminum Alloy

Function: Designed as an accessory for microscopes to modify magnification and enhance imaging capabilities.

Magnification: 0.5x, reducing the standard magnification to provide a wider field of view.

Working Distance: 165 mm

Compatibility: Fits most standard stereo microscopes, enhancing flexibility in magnification options.

Durability: Robust construction with an aluminum alloy body, ensuring long-lasting performance.

Weight: Lightweight design for ease of use and minimal impact on the microscope's balance.

Package List: Includes 1 unit 0.5X Auxiliary Objective Lens in the package.

Application: Ideal for applications requiring lower magnification for broader viewing areas, such as inspection, assembly, or dissection tasks.

