



PZ-3
with BULL LENS

Dissecting Microscope

with BULL LENS

- Senior Dissecting Microscope Has Heavy Round Base With A Precisely Designed Body. Sensitive Focusing Is Done By Rack-pinion Movement. A Joint Revolving Arm Is Provided For Moving Magnifying Lens Overfull Area Of The Stage. It Is Easily Replaceable By A Special Bulls Eye Lens (83mm Dia.) To Increase The Field Of Vision.

- A Plano-concave Mirror Reflects Sufficient Day Light And Artificial Light On The Object. Detachable Hand Rests On Both Sides Of Stage Provided For A Study Support During Dissection. Provided With 10x & 20x Eye Pieces And Special Bulls Eye Lens (83mm Dia.).

Packaging: Inner Thermocol/Styrofoam Packing And Outer Cardboard Box.

