



QHD-SZ

Stereo Zoom Microscope

With SCREEN and Inbuilt Camera

Display: IPS 10-inch Touch Screen, Adjustable From 5 To 15 Degrees.

Camera: Built-in HDR image, OSD Manual, 1/3" CMOS Sensor, 4.0 Mega Pixels, 16GB Memory.

Zoom Lens: 0.7x To 5x Magnification Range.

Output Options: HDMI/USB2.0/Mini USB.

Storage Support: SD Card.

Peripheral Compatibility: Wireless Keyboard And Mouse Support.

Main Objective: 1x/0.025 Main Objective With A Working Distance (WD) Of 105mm.

Illumination:

- LED Ring Lamp With 56 Bulbs, Adjustable Intensity, Colour Temperature 5000-5000K.
- Two-circle LED ring lamp, separately adjustable intensity, colour Temperature 5000-5000K.

Base: Big plain base, lifting range 225mm.

Power: Universal Power Supply 110v-240v Integrated With Voltage Flux Stabilizer.

Packaging: Inner Thermocol / Styrofoam Packing And Outer Cardboard Box.

