



SM-50

Binocular Stereo Zoom Microscope

Optical System: Greenough Optical System.

Stand: Ergonomically Designed Metal Pillar Stand.

Viewing Head: Binocular Head 45° Inclined 360° Rotatable, IPD 48-75mm With Diopter Adjustment.

Eye Guards: Adjustable Rubber Eye Guards For Comfortable Viewing.

Eyepieces: Wide Field 10x, F.O.V 18mm With ± 5 Dioptic Adjustment, Antifungal & Anti-reflection Coated.

Objectives: 2X/0.05, 4x/0.10 Revolving On Turret.

Magnification: 20x To 40x.

Working Distance: 60mm.

Stand: Dual light stand.

Focusing: Coaxial Coarse And Fine Focusing. With Fine Adjustment Scale Value 0.002mm.

Illumination: 3.5W LED Light (Up To 100,000 Working hours) with Intensity Controller.

Power: Universal Power Supply 110v-240v Integrated With Voltage Flux Stabilizer.

Packaging: Inner Thermocol / Styrofoam Packing And Outer Cardboard Box.

