

Trinocular Stereo Zoom Microscope With Beam Splitter

Optical System: Greenough Optical System.

Stand: Ergonomically Designed Metal Pillar Stand.

Viewing Head: Binocular Head 20° Inclined 360° Rotatable, IPD 48-75mm With Diopter Adjustment and Beam Splitter.

Eye Guards: Adjustable Rubber Eye Guards For Comfortable Viewing.

Eyepieces: Wide Field 10x, F.O.V 22mm With ± 5 Diopteric Adjustment, Antifungal & Anti-reflection Coated.

Magnification Range: 8x To 80x Continuous.

Zoom Ratio: 1 : 10

Working Distance: 110mm.

Lenses: Anti-fungus Treated And Anti-reflection Coated.

Focusing: Coaxial Coarse And Fine Focusing. With Fine Adjustment Scale Value 0.002mm.

Illumination: Dual 3W Incident & Transmitted LED Light (Up To 100,000 Working hours) Combined Illumination Control Buttons With Rheostat Dial For Light Intensity Control.

Power: Universal Power Supply 110v-240v Integrated With Voltage Flux Stabilizer.

Packaging: Inner Thermocol/Styrofoam Packing And Outer Cardboard Box.

