



SZB-45E

Trinocular Stereo Zoom Microscope

Optical System: Greenough Optical System.

Stand: Single Mould Sturdy Aluminum Die-cast Stand (Anti-Rust) Body.

Viewing Head: Trinocular Head 45° Inclined 360° Rotatable, IPD 48-75mm With Diopter Adjustment.

Eye Guards: Adjustable Rubber Eye Guards For Comfortable Viewing.

Eyepieces: Wide Field 10x, F.O.V 20mm With +/-5 Diopter Adjustment, Antifungal & Anti-reflection Coated.

Magnification Range: 7x To 45x Continuous.

Zoom Ratio: 1 : 6.4

Working Distance: 110mm.

Illumination: Dual 3W Incident & Transmitted LED Light (Up To 100,000 Working hours) Combined Illumination Control Buttons With Rheostat Dial For Light Intensity Control.

Focusing: Coaxial Coarse And Fine Focusing. With Fine Adjustment Scale Value 0.002mm.

Power: Universal Power Supply 110v-240v Integrated With Voltage Flux Stabilizer.

Packaging: Inner Thermocol/Styrofoam Packing And Outer Cardboard Box.

